


Products Upgrade Announcement

All PakGent production lines have fully been upgraded especially on packaging, COA, and label etc.



Part 1. Pipette tips rack upgrade




Rack weight increased by **33%**




Jet coding information
displayed on the rack



Lids can be opened
180 degrees



Buckle thickness
increased by **32%**



Rack pack thickness



▲ After

▲ Before

Label information



▲ After

▲ Before

Buckle thickness



▲ After

▲ Before

Lid opening Angle



▲ After

▲ Before



UFPT-F-10	UFPT-F-10E	UFPT-F-10ET
UFPT-10R	UFPT-10ER	UFPT-10ETR
UFPT-F-10-L	UFPT-F-10E-L	UFPT-F-10ET-L
UFPT-10R-L	UFPT-10ER-L	UFPT-10ETR-L

UFPT-F-20	UFPT-F-50	UFPT-F-100
UFPT-F-200	UFPT-F-200S	UFPT-200R
UFPT-200SR	RT20R	RT20RF
RT200R	RT200RF	UFPT-F-20-L
UFPT-F-50-L	UFPT-F-100-L	UFPT-F-200-L
UFPT-F-200S-L	UFPT-200R-L	UFPT-200SR-L
RT20R-L	RT20RF-L	RT200R-L
RT200RF-L		

UFPT-F-1000E	UFPT-F-1000M
UFPT-1000ER	UFPT-1000MR
RT1000RF	RT1200RF
RT1000R	RT1200R
UFPT-F-1000E-L	UFPT-F-1000M-L
UFPT-1000ER	UFPT-1000MR-L
RT1000RF-L	RT1200RF-L
RT1000R-L	RT1200R-L

Part 2. COA upgrade (All products)

PakGent®
 @ www.pakgentbio.com
 @ info@pakgentbio.com

Certificate of Analysis

This is to certify that units of the lot number below were tested and found to comply with the specifications of the grade listed. Products are tested to meet the analytical requirements of the noted grade. The following information is the actual analytical results obtained, for reference.

Part Number:	UFPT-200R-L		
Lot.No.:	202207300101		
Specification:	Volume:	<input type="checkbox"/> 0.1-10ul <input type="checkbox"/> 0.1-20ul <input type="checkbox"/> 1-20ul <input type="checkbox"/> 1-50ul <input type="checkbox"/> 1-100ul <input checked="" type="checkbox"/> 1-200ul <input type="checkbox"/> 1-300ul <input type="checkbox"/> 100-1000ul <input type="checkbox"/> 100-1200ul <input type="checkbox"/> 0.5-5ml <input type="checkbox"/> 1-10ml	
	Filter:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Package: <input type="checkbox"/> bag <input checked="" type="checkbox"/> rack <input type="checkbox"/> refill	
Product Name:	Pipette tips(Low Retention)		
Package:	96pcs/rack, 50racks/case		
Shelf Life:	3 years		
Manufacturer:	PakGent Bioscience (Suzhou) Co., Ltd.		
Inspection Date:	2th,February,2024	Report Date	2th,February,2024
Inspection items and result:			
Cosmetic:	No black spots, no cuts,clear and transparent.		
Package:	Accurate quantity and label, inner package without damage		
Dimension:	OK		
Function test:	Accuracy≥98%		
Endotoxin:	<0.03EU/ml		
RNase/DNase:	DNase≤4x10 ⁻⁶ Kunitz units/ul, RNase≤10 ⁻⁶ Kunitz units/ul		
ATP:	2×10 ⁻¹² mg/ul		
PCR Inhibitors:	Not detectable		
Human DNA:	<30pg		
Sterile:	Irradiation,SAL10 ⁻⁶		
Conclusion: Comply with Specification			

This signature indicates that the above test has met all quality specification and has been reviewed by a quality representative.

Title: QA Manager
Signature: Alvin Cao



Visible clarity in certification

Accuracy:	Accuracy ≥98%
Endotoxin:	<0.03EU/ml
RNase:	DNase≤4x10 ⁻⁶ Kunitz units/ul
ATP:	2×10 ⁻¹² mg/ul
PCR Inhibitors:	Not detectable
Human DNA:	<30pg
Sterile:	Irradiation,SAL10 ⁻⁶

Add indicator parameters

PakGent®
 @ www.pakgentbio.com
 @ info@pakgentbio.com

Certificate of Analysis

This is to certify that units of the lot number below were tested and found to comply with the specifications of the grade listed. Products are tested to meet the analytical requirements of the noted grade. The following information is the actual analytical results obtained, for reference.

Part Number:	UFPT-200R-L		
Lot.No.:	202207300101		
Specification:	Volume:	<input type="checkbox"/> 0.1-10ul <input type="checkbox"/> 0.1-20ul <input type="checkbox"/> 1-20ul <input type="checkbox"/> 1-50ul <input type="checkbox"/> 1-100ul <input checked="" type="checkbox"/> 1-200ul <input type="checkbox"/> 1-300ul <input type="checkbox"/> 100-1000ul <input type="checkbox"/> 100-1200ul <input type="checkbox"/> 0.5-5ml <input type="checkbox"/> 1-10ml	
	Filter:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No Package: <input type="checkbox"/> bag <input checked="" type="checkbox"/> rack <input type="checkbox"/> refill	
Product Name:	Pipette tips(Low Retention)		
Package:	96pcs/rack, 50racks/case		
Shelf Life:	3 years		
Manufacturer:	PakGent Bioscience (Suzhou) Co., Ltd.		
Inspection Date:	2th,February,2024	Report Date	2th,February,2024
Inspection items and result:			
Cosmetic:	No black spots, no cuts,clear and transparent.		
Package:	Accurate quantity and label, inner package without damage		
Dimension:	OK		
Function test:	Accuracy≥98%		
Endotoxin:	<0.03EU/ml		
RNase/DNase:	DNase≤10 ⁻⁶ Kunitz units/ul, RNase≤10 ⁻⁶ Kunitz units/ul		
ATP:	2×10 ⁻¹² mg/ul		
PCR Inhibitors:	Not detectable		
Human DNA:	<30pg		
Sterile:	Irradiation,SAL10 ⁻⁶		
Conclusion: Comply with Specification			

This signature indicates that the above test has met all quality specification and has been reviewed by a quality representative.

Title: QA Manager
Signature: Alvin Cao

▲ After

PakGent

Certificate Of Analysis

Product Name	Filter Tips
Lot No.	202207080101
Specification:	1-200ul (UFPT-F-200-L)
Packing	96/rack, 50racks/case
Shelf life	3 years
Origin	PakGent Bioscience (suzhou) Co.,Ltd.
Inspection Date	8th, July, 2022
Report Date	8th, July, 2022
Appearance:	No black spots, no cuts, no cuts, clear as a bell, clear and transparent
Package:	Accurate quantity and label, inner packing without damage
Performance test:	Accuracy ≥99%
Endotoxin (pyrogen) Testing:	
The product has been tested and has met the criteria established in the Guideline on Validation of the Limulus Amebocyte Lysate Test as an End-Product Endotoxin Test for Human and Animal Parenteral Drugs, Biological Products and Medical Devices (FDA, December, 1987). No pyrogens were detected.	
RNase/DNase Testing:	
This product has been tested and is free of any detectable RNase/DNase contamination.	
Test/Procedure	Result
Visual Attributes	Pass
Dimensional	Pass
Performance	Pass
Sterile	Pass
Test Result: Up to Specification	
This signature indicates that the above test has met all quality specification and has been reviewed by a quality representative	
Title: QA Director	
Signature: Emily	

▲ Before

Part 3. The outer case label upgrade (All products)

LTS Tip Visible clarity in Certification

STERILE Sterility Grade SAL 10 Visible clarity in sterilization information

Expiry Date Visible clarity in expiration date

PakGent Bioscience Available contact information

RT200RF-L LTS Tip

200ul
LTS tip
with filter, low retention
DNase/RNase free
96/rack, 50racks/case


STERILE Sterility Grade SAL 10
If packaging sealed and undamaged

Expiry Date 15/09/2026

PakGent Bioscience (SuZhou) Co.,Ltd
+86 (0) 512 6670 0017
info@pakgentbio.com
www.pakgentbio.com



▲ After



▲ Before

At present, we are working to switch to the new production lines. In this process, both new and previous products will exist simultaneously for a short period of time. We would like to express our sincere apologies for any inconvenience caused.

Thank you again for your support of our product, and we look forward to your valuable feedback. Wish you all the best!