

Technical Datasheet

Article - No:	TDS-PCRS-10F	Revision: A0																																				
Remarks	single use only																																					
Specification																																						
Product name	PCRS-10F, 0.1ml PCR 8-strip qPCR tube																																					
Production conditions	Class 100,000 clean room production																																					
	No release agent is used																																					
Classification	CE/IVD																																					
Dimensions	<table border="1" style="margin-top: 10px; width: 100%; text-align: center;"> <tr> <td>A1</td><td>B1</td><td>C1</td><td>D1</td><td>E1</td><td>F1</td><td>G1</td> </tr> <tr> <td>5.56</td><td>76.91</td><td>62.93</td><td>7.38</td><td>0.65</td><td>5.43</td><td>5.29</td> </tr> <tr> <td>A2</td><td>B2</td><td>C2</td><td>D2</td><td>E2</td><td>F2</td><td>G2</td> </tr> <tr> <td>5.52</td><td>80</td><td>62.6</td><td>7.3</td><td>0.74</td><td>5.43</td><td>5.58</td> </tr> <tr> <td>1</td><td>2</td><td>3</td><td>5</td> </tr> <tr> <td>80</td><td>64</td><td>15</td><td>1.3</td> </tr> </table>		A1	B1	C1	D1	E1	F1	G1	5.56	76.91	62.93	7.38	0.65	5.43	5.29	A2	B2	C2	D2	E2	F2	G2	5.52	80	62.6	7.3	0.74	5.43	5.58	1	2	3	5	80	64	15	1.3
A1	B1	C1	D1	E1	F1	G1																																
5.56	76.91	62.93	7.38	0.65	5.43	5.29																																
A2	B2	C2	D2	E2	F2	G2																																
5.52	80	62.6	7.3	0.74	5.43	5.58																																
1	2	3	5																																			
80	64	15	1.3																																			
Colour	Clear																																					
Nominal Capacity	0.1ml																																					
Material Characteristics	polypropylene (PP), premium surface																																					
Product Characteristics	biologically inert																																					
	no heavy metals according to EC/94/62																																					
Purity Grade	free of RNase and DNase																																					
	free of pyrogen																																					
	free of PCR inhibitors																																					
	free of human DNA																																					
	free of endotoxins																																					
Shelf - Life	60 months at controlled storage conditions																																					
Packaging	pcs per box: 125																																					
	box per case: 10																																					
Storage	dry condition: 40%-60% rel. humidity																																					
Approved by PakGent quality department																																						